FILLMED
LABORATOIRES

NCTF[®]
135 HA

USE

SYMBOLS ON THE PACKAGING

CE
0297

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CONDITIONS DE STOCKAGE/CONSERVATION

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NCTF® 135 HA
Highly concentrated hyaluronic acid
+ polyrevitalising solution

INSTRUCTIONS FOR USE

DESCRIPTION

Kit of 5 vials

5 x 3 mL vials
5 x 3 mL syringes
5 x 18G withdrawal needles
5 x 30G injection needles
5 x 32G injection needles

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5 nucleic acids

Deoxyadenosine
Deoxycytidine
Deoxyguanosine
Deoxythymidine
5-Methyl-2'-deoxycytidine

24 amino acids

α-Aminobutyric acid
Alanine
Arginine
Asparagine
Aspartic acid
Cystine
Glutamine
Glutamic acid
Glycine
Histidine
Hydroxyproline
Isoleucine
Leucine
Lysine
Methionine
Ornithine
Phenylalanine
Proline
Serine
Taurine
Threonine
Tryptophane
Tyrosine
Valine

6 coenzymes

TPP (Coccarboxylase)
CoA (Coenzyme A)
FAD (Flavine adenine dinucleotide)
NAD (Nicotinamide adenine dinucleotide)
NADP (Nicotinamide adenine dinucleotide phosphate)
UTP (Uridine triphosphate)

Others

Glutathione
Polysorbate 80
Glucuronic acid
Glucuronic acid lactone
Glucosamine
Dextrose anhydrous

INDICATIONS

NCTF® 135 HA is an anti-ageing injectable product indicated for intense revitalization and hydration of tired or dull skin, filling of superficial wrinkles and redensification of mature or slack skin. NCTF® 135 HA is indicated to be injected in the epidermis or dermis in the following areas:

- face (including periorbital area)
- neck
- décollete area
- back of the hands

- internal face of the arms

CONTRAINDICATIONS

- Do not inject into blood vessels. Do not use in:
- patients with a skin alteration or disease of any kind in the treated area;
 - patients with a history of auto-immune disease or under immunotherapy;
 - patients with hypersensitivity or a known allergy to any of the constituents;
 - pregnant or breast-feeding women;
 - children.

PRECAUTIONS FOR USE

There are no clinical data available in terms of tolerance with regard to the injection of NCTF® 135 HA in an area that has already been treated with another "permanent" filling product. Precaution should be taken in such cases, although NCTF® 135 HA is not injected at the same level. Patients receiving an anticoagulant treatment or a platelet aggregation inhibitor (for example, aspirin) must be informed of the increased risk of bruising and mild bleeding at the injection site. The patient should be given the recommendation not to use makeup for 12 hours following the injection and to avoid prolonged exposure

to the sun, UV rays, freezing temperatures as well as use of saunas or steam rooms for at least 48 hours after the injection. Training in injection techniques is recommended for the use of NCTF® 135 HA. NCTF® 135 HA must not be injected at the same time as a medium or deep chemical peeling or dermabrasion. A local anaesthetic (patch or cream) can be applied prior to the injection. Please note that anaesthetic products can cause local redness or hypersensitivity. Do not inject NCTF® 135 HA into beauty spots (naevi) and in active acne lesions.

INCOMPATIBILITIES

There is a known incompatibility between hyaluronic acid and quaternary ammonium salts such as benzalkonium chloride. NCTF® 135 HA should therefore never come into contact with such products or with medical/surgical equipment treated with this type of product.

UNDESIRABLE EFFECTS

The practitioner must inform the patient that there are possible undesirable effects related to the implantation of this device, such as: redness or mild local inflammation that usually disappears within 24 to 48 hours; mild oedema and small haematomas may occur, but disappear within 48 hours; minimal bleeding at the injection site that stops rapidly and spontaneously after the injection; transient pain at the injection site. The patient must notify the practitioner as soon as possible of any inflammatory reactions lasting more than one week or the appearance of any other adverse effect. The practitioner must treat the adverse effect with suitable treatment. Any other adverse effect related to the injection of NCTF® 135 HA must be reported to the distributor and/or manufacturer.

INSTRUCTIONS FOR USE

Before starting treatment, the patient must be informed of the indications of the device, its contraindications, its incompatibilities and its possible undesirable effects. All asepsis and hygiene precautions must be taken to prepare the site, the material and the patient prior to the injection. NCTF® 135 HA injection can be performed manually or in association with a device like an injection gun.

Kit of 5 vials

- 1/ Connect a 18G needle to a 3 mL syringe (provided in the kit) and aseptically withdraw NCTF® 135 HA solution to be injected.
- 2/ Disconnect the 18G needle and connect the suitable sterile needle for injection provided in the kit (30G or 32G) to the syringe filled with NCTF® 135 HA.
- 3/ Remove the cap of the needle and inject NCTF® 135 HA.

Box of 10 vials

NCTF® 135 HA must be injected with sterile, single use material (needles for micro-injection), which is not provided.

Data relating to the injection (patient details and product batch number) must be kept by the practitioner.

PROTOCOL

The recommended dose is 0.05 mL per injection point, spaced around 1 – 1.5 cm. The recommended protocol consists in 3 treatments, at intervals of 3-4 weeks. Additional treatments can be given every 3 or 6 months if the practitioner deems it necessary with regard to skin quality.

WARNINGS

- Check the expiry date on the label.
- Do not reuse.
- Do not resterilize.
- After use, the needles must be thrown away in a suitable safe container.
- Reuse of the product involves risks (e.g. cross contamination) for the patient.
- Do not use if the packaging is open or damaged.

**STORAGE/CONSERVATION CONDITIONS**

Store NCTF® 135 HA between 2 and 25°C (36-77°F) in a dry place protected from light and moisture.
Handle with caution.

**QUICK AUTHENTICITY CHECK**

Thank you for buying a FILL-MED product.

You can check this product's authenticity using our special protection: the Bubble Tag™ ③ on the ProofTag placed on the outer packaging.

You have two options:

- option A: scan the QR code ①,
- option B: enter the authentication number ② on www.proof.fillmed.com.

PROOFTAG

- ① QR code to scan with your smartphone.
- ② Authentication number.
- ③ Bubble Tag™ to compare with the information available on the www.proof.fillmed.com website.

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