FILLMED
LABORATOIRES

NCTF[®]
135 HA

USE

NCTF®
135 HA

SYMBOLS ON THE PACKAGING

0477

THE NEEDLES AND SYRINGES CONTAINED
IN THE PACK ARE CE MARKED UNDER
THE RESPONSIBILITY OF THE MANUFACTURER
(SEE LABELS)

MANUFACTURER

DATE
OF MANUFACTURE

DO NOT RESTERILIZE

DO NOT RE-USE

STORAGE
TEMPERATURE

CAUTION

PROTECT FROM
SUNLIGHTREAD THE NOTICE/
THE INSTRUCTIONS
BEFORE USE

SINGLE STERILE BARRIER SYSTEM

THE PACKAGING MUST
REMAIN DRY

EXPIRY DATE

STERILE

A

STERILIZED WITH ASEPTIZED
TECHNIQUESDO NOT USE
IF THE PACKAGING
IS DAMAGED

LOT

BATCH NUMBER

QUICK AUTHENTICITY CHECK

③

①

②

Thank you for buying a FILL-MED product.
You can check this product's authenticity using our special protection:
the Bubble Tag™ ③ on the Proof tag placed on the outer packaging.

You have two options:

- option A: scan the QR code ①,
- option B: enter the authentication number ② on www.proof.fillmed.com.

PROOFTAG

- ① QR code to scan with your smartphone.
- ② Authentication number.
- ③ Bubble Tag™ to compare with the information available on the www.proof.fillmed.com.

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CONDITIONS DE STOCKAGE/CONSERVATION

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HEALTHCARE PROFESSIONALS LEAFLET
Only to be administered by appropriately trained healthcare professionals who are qualified or accredited in accordance with national law.

NCTF® 135 HA

Highly concentrated hyaluronic acid + polyrevitalising solution

DESCRIPTION

Kit of 5 vials

5 x 3 mL vials
5 x 3 mL syringes
5 x 18G withdrawal needles
5 x 30G injection needles
5 x 32G injection needles

COMPOSITION

- **Non-crosslinked sodium hyaluronate / Hyaluronic acid** (5 mg/mL) (constituent responsible for the intended action) – Average molecular weight: 0.7 to 1.3 x 10⁶ Dalton
Function: filling effect and hydration
- **Buffer solution** (final pH 6.50-7.80)

Function: adjustment of pH

- **Polyrevitalising solution** (9 mg/mL)

Function: provision of an optimal physiological environment for skin cells

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6 minerals: calcium chloride, potassium chloride, magnesium sulfate, sodium acetate, sodium chloride, sodium dihydrogenophosphate

5 nucleosides: deoxyadenosine, deoxycytidine, deoxyguanosine, deoxythymidine, 5-methyl-2'-deoxycytidine

24 amino acids: α-aminobutyric acid, alanine, arginine, asparagine, aspartic acid, cystine, glutamine, glutamic acid, glycine, histidine, hydroxyproline, isoleucine, leucine, lysine, methionine, ornithine, phenylalanine, proline, serine, taurine, threonine, tryptophane, tyrosine, valine

6 coenzymes: TPP (cocarboxylase), CoA (coenzyme A), FAD (flavine adenine dinucleotide), NAD (nicotinamide adenine dinucleotide), NADP (nicotinamide adenine dinucleotide phosphate), UTP (uridine

triphosphate)

Others: glutathione, polysorbate 80, glucuronic acid, glucuronic acid lactone, glucosamine, dextrose anhydrous

INTENDED USE /INDICATIONS

NCTF® 135 HA products are not for medical use but for aesthetic purposes. NCTF® 135 HA is an anti-ageing injectable product indicated for filling of superficial wrinkles for intense revitalization and hydration of tired or dull skin, and to improve skin homogeneity and radiance. It also helps skin redensification, thickness and elasticity of mature or slack skin. NCTF® 135 HA is indicated to be injected in the dermis in the following areas:

- face (including periorbital area)
- neck
- décolleté area.

CONTRAINDICATIONS

Do not inject into blood vessels. Do not use in:

- patients with a skin alteration or disease of any kind in the treated area;
- patients with a history of auto-immune disease or under immunotherapy;
- patients with hypersensitivity or a known allergy to any of the constituents;
- patients who are less than 18 years old
- pregnant or breast-feeding women.

PRECAUTIONS FOR USE

Do not use the product if the product aspect

(color, transparency...) changed. There are no clinical data available in terms of tolerance with regard to the injection of NCTF® 135 HA in an area that has already been treated with another "permanent" (no resorbable) filling product. Precaution should be taken in such cases, although NCTF® 135 HA is not injected at the same level. Patients receiving an anticoagulant treatment or a platelet aggregation inhibitor (for example, aspirin) must be informed of the increased risk of bruising and mild bleeding at the injection site. The patient should be given the recommendation not to use makeup for 12 hours following the injection and to avoid prolonged exposure to the sun, UV radiation, hot and cold temperatures as well as use of saunas or steam rooms for at least 48 hours after the injection. Training in injection techniques is recommended for the use of NCTF® 135 HA by healthcare professionals. NCTF® 135 HA must not be injected at the same time as a medium or deep chemical peeling or dermabrasion. A local anaesthetic (patch or cream) can be applied prior to the injection. Please note that anaesthetic products can cause local redness or hypersensitivity. Do not inject NCTF® 135 HA into beauty spots (naevi) and in active acne lesions.

INCOMPATIBILITIES

Hyaluronic acid is incompatible with quaternary ammonium compounds, such as

benzalkonium chloride solutions. This is why NCTF® 135 HA must never come into contact with medical and surgical instruments treated with this type of product.

ADVERSE EFFECTS

The healthcare professional must inform the patient that there are possible undesirable effects related to the implantation of this device, such as: allergic systemic reaction (angioedema, allergic reaction, itching sensation, rash, chills, fever), local inflammatory reaction (erythema, swelling, burning sensation, nodule, granuloma, pain, skin inflammation), local allergic reaction (pruritus) and infection (abscess, herpes). The healthcare professional must ensure that the patient notify them of any undesirable effects.

DECLARATION OF ADVERSE EFFECTS

The patient must inform the healthcare professionals as soon as possible if any side effect (cited above or other) occurs or if any side effect (cited above or other) persists for longer than one week. The implant must be checked and treated appropriately. The healthcare professionals will treat these adverse effects appropriately (see next section). The distributor and/or manufacturer must be alerted of any other adverse effect related to the injection of NCTF® 135 HA. You can report adverse reactions at the manufacturer directly by contacting: vigilance@fillmed.com or on

the Laboratoires FILL-MED's website: <https://fillmed.com>. Any serious incident must be reported to the authorities of the member state in which the user and/or patient is established. The summary of safety and clinical performance (SSCP) is available on Laboratoires FILL-MED's website.

TREATMENTS OF ADVERSE EFFECTS

- Inflammatory reaction (swelling, erythema, rash, redness, nodule, granuloma, skin inflammation...) can be treated by: cold compress for 48 hours, anti-inflammatory topical creams, sun creams, antibiotics, and corticosteroids if necessary.
 - Skin irregularities (nodule, papula) can be treated with hyaluronidase or massage of the involved zone.
 - Allergic reaction (itching sensation, burning sensation) can be treated by: antihistamines and corticosteroids and hospitalization in case of angioedema/quincke edema.
 - Infection (abscess) can be treated by antibiotics and surgical intervention, for cold sores (herpes) antiviral medication.
 - Systemic adverse events (chills, fever...) can be treated by: medical intervention by health professional.
 - Hospitalization (in severe cases).
- The prolonged use of any medication in treatment of adverse events has to be carefully assessed, since this may carry a

risk for the patient.

INSTRUCTIONS FOR USE

NCTF® 135 HA is intended to be injected into the dermis by a healthcare professional trained in injection techniques. The healthcare professional can use injection techniques adapted to the area being treated with needle, micro needle, or canula, such as: multiple intradermal injection/ point per point, retracing/ withdrawal (simple or fanning), nappage. Before starting treatment, the patient must be informed of the indications of the device, its contraindications, its incompatibilities and its possible undesirable effects. All asepsis and hygiene precautions must be taken to prepare the site, the material and the patient prior to the injection. NCTF® 135 HA injection can be performed manually or in association with sterile and single-use syringes and needles in accordance with their intended use (intra-dermal injections). After the treatment, the healthcare professional must provide the patient with a product leaflet, implant card on which the traceability label must be stuck. Data relating to the injection (patient details and product batch number) must also be kept by the healthcare professional.

Kit of 5 vials

1/ Connect a 18G needle to a 3 mL syringe (provided in the kit) and aseptically

withdraw NCTF® 135 HA solution to be injected.

2/ Disconnect the 18G needle and connect the suitable sterile needle for injection provided in the kit (30G or 32G) to the syringe filled with NCTF® 135 HA.

3/ Remove the cap of the needle and inject NCTF® 135 HA.

PROTOCOL

The recommended dose is 0.05 mL per injection point, spaced around 1 – 1.5 cm. The recommended protocol consists in 3 treatments, at intervals of 3-4 weeks. Additional treatments can be given every 3 or 6 months if the healthcare professional deems it necessary regarding the skin quality. The maximum injected quantity is 6 mL for the face, 3 mL for the neck and 3 mL for the décolleté area per session. Free hyaluronic acid resorption time is expected to be between 24 and 48 hours.

WARNINGS

- Check the expiry date on the label.
- Check the integrity of the sterile barrier immediately prior to using the product.
- Do not reuse the vial & the components on other patients.
- Do not resterilize.
- After use, the needles must be thrown away in a suitable safe container.
- Reuse of the product involves risks (e.g. cross contamination) for the patient.

- Do not use if the packaging is open or damaged.
- Do not inject excessive volume, it can potentially lead to side effects.

STORAGE/CONSERVATION

CONDITIONS

Store NCTF® 135 HA between 5°C and 25°C (41-77°F) in a dry place protected from light and moisture. Handle with caution.

STANDARDS

The product is certified CE and compliant with harmonised standards and common specifications (EU) 2022/2346 applicable to NCTF® 135 HA referenced in SSCP-NCTF-NCL.

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